

Leckey Pronestander



- Assists in standing and weight bearing development
- Adjustable in 4 forward prone positions up to 40 degrees
- Highly adjustable to suit individual users
- Includes chest and knee supports, hip belt, angle adjustable tray and activity bowl
- Mobility base has 4 locking castors for stability
- User Standing Height - PED479250 = 800 - 1150 mm, PED479260 = 1100 - 1450 mm, PED479270 = 1400 - 1850 mm
- Chest Pad Height - PED479250 = 470 - 680 mm, PED479260 = 680 - 975 mm, PED479270 = 975 - 1350 mm



ndis

Item Code:

05_053603010_0103_1_2

Click [here](#) for more information

Warranty Information

- 2 Years

LECKEY®

Code	Display Name	Safe Working Load	Unit Weight
PED479250	Leckey Pronestander - Size 1, Ages 1.5-5 years	60 kg	34 kg
PED479260	Leckey Pronestander - Size 2, Ages 4-10 years	80 kg	38 kg
PED479270	Leckey Pronestander - Size 3, Ages 9-18 years	90 kg	49 kg