

KeraSit2Sit Accessories



- **Backstrap** – Specially designed for the Kera sit2sit, a One Size Fits Most backstrap will fit most users comfortably. The Kera sit2sit comes with one Backstrap included.
- **Thighstrap** – Using the Thighstrap improve the comfort of the user by allowing reduced tightness from the Backstrap during transfers. It can also reduce any weight bearing through the user's feet. Easily attached to any Kera sit2sit.
- **Footbox** - Provide a means of ensuring the user's feet are held securely in place during a transfer. Come as a pair and include a pair of Footbox Straps to secure the behind the heel or over the top of the foot. Easily fitted to the baseframe of the Kera sit2sit.
- **Footbox Straps** - Additional pair of straps for use with Kera Footboxes to further secure the user's feet.
- **Riser Blocks** - Foam blocks are easily fitted into the Footboxes and raise the platform when the user's feet don't reach the Footplates. Two heights available: 60mm or 110mm.
- **Knee Foam Pack** – Foam inserts to extend the knee surface to contact the user's knee if needed. Foam packers are designed to fit behind the Knee Support Covers, four 11mm thick inserts allow 11mm-22mm extension on both knees.
- **Petite Chestpad** – Reduced size Chestpad for smaller users or children. Can be swapped out for the standard size Chestpad.



Contract approved

Support @ Home

Function: Changing body position

Category: Prescribed

Click [here](#) for more information

Code	Display Name
LSS398300	Backstrap - Kera Sit2Sit
LSS398310	Thighstrap - Kera Sit2Sit
LSS398320	Footbox - Kera Sit2Sit
LSS398330	Footbox Straps - Kera Sit2Sit
LSS398340	Knee Foam Pack - Kera Sit2Sit
LSS398350	Backstrap extenders - (pair) - Kera Sit2Sit
LSS398360	Hook Lowering Bracket - Kera Sit2Sit
LSS398370	Chest Pad Assembly - Standard - Configurable - (includes cradle & slider bar)
LSS398380	Chest Pad Assembly - Short - Configurable - (includes cradle & slider bar)