

The Medication Carts



- Adjustable push handle, quick and easy operation to enable users to set the push handle at the correct ergonomic height
- Lightweight composite aluminium/plastic panels, reduces cart noise while in motion and panels can be replaced economically if damaged
- Recessed storage compartment
- Recessed ABS formed plastic work surface
- Drawers fitted with soft close runners and damper - prevents drawers from opening while in motion
- Anti-static total lock castors
- Raised drawer handles designed for grip and easy cleaning
- Corner buffers
- Epoxy powder coated frame



Warranty Information

- 5 Years Frame

Code	Display Name	Total Length	Total Width	Total Height
WAT901140	Medication Cart - 5 Drawer + Dividers (12)	-	-	-
WAT901130	Medication Cart - 4 Drawer + Dividers (12)	-	-	-
WAT901120	Medication Cart - 3 Drawer (Webster)	-	-	-
WAT901150	Med/Theatre Wide Cart - 4 Drawer + Dividers (12 & 14)	570 mm	870 mm	1065 mm

Code	Display Name	Total Length	Total Width	Total Height
WAT901160	Med/Theatre Wide Cart - 5 Drawer + Dividers (12 & 14)	570 mm	870 mm	1065 mm