

## Over Toilet Aid Maxi - 600mm width (300kg)



- Extra wide and reinforced mild steel tubing frame accomodating maximum user weights up to 300 kg
- Fixed height frame for further strength and stability
- Removable lift back seat for easy cleaning, comes standard with splash guard
- Fitted with non-slip rubber tips
- Zinc treated anti corrosive frame



\*\*\* Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

Code	Display Name
H_BTT146500	Over Toilet Aid Maxi - 600mm width (300kg)