

McDavid Neoprene Ligament Support



- Delivering firm compression and soft tissue support.
- Made of medical grade lightweight neoprene that retains body warmth.
- 4 elastic straps that deliver additional support to ligaments.
- Horse-shoe buttress to stabilise the knee cap.
- contoured design for a better fit with two medial and two lateral medical steel stays.
- Top and bottom anchor straps minimise slippage.

✓ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Colour	Knee Circumference
FBS284860	McDavid Neoprene Ligament Support - Black - Sml	Black	330 - 356 mm
FBS284870	McDavid Neoprene Ligament Support - Black - Med	Black	356 - 381 mm
FBS284880	McDavid Neoprene Ligament Support- Black - Lge	Black	381 - 406 mm
FBS284881	McDavid Neoprene Ligament Support - Black - XL	Black	406 - 445 mm
FBS284882	McDavid Neoprene Ligament Support - Black - XXL	Black	445 - 483 mm