


Aspire Active Air Bari 10 Bariatric Mattress Replacement



- Integrated Turning Therapy - Allows patient rotation up to 40° for improved pressure care and to assist staff with patient handling requirements
- Rapid Pump Inflation - Maximises therapy response times and minimises staff resources by fully inflating in only 2 minutes
- Hospital Grade Cover – Thicker weight polyurethane fabric with increased bleach resistance to handle modern hospital infection control standards. Suitable for Sodium Hydrochloride solutions
- Adaptable Pump - Pump can be switched between modes for typical users (up to 160 kg) and bariatric users (up to 450 kg) ensuring a single asset can meet the broadest spectrum of patient needs

 Contract approved

Compliance

IEC 60601-1-2

Warranty Information



- Cover: 2 Years
- Air Mattress: 2 Years
- Pump: 2 Years

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Safe Working Load	Pressure Care Rating
PTM639500	Active Air BARI 10 - Single	12019	1	2000 mm	900 mm	270 mm	30 kg - 450 kg	Very High Risk
PTM639510	Active Air BARI 10 - King Single	12679	1	2000 mm	1070 mm	270 mm	30 kg - 450 kg	Very High Risk
PTM639520	Active Air BARI 10 - Super King Single	13339	1	2000 mm	1230 mm	270 mm	30 kg - 450 kg	Very High Risk