

Aspire Active Air Bari 10 Bariatric Mattress Replacement




- Integrated Turning Therapy Allows patient rotation up to 40° for improved pressure care and to assist staff with patient handling requirements
- Rapid Pump Inflation Maximises therapy response times and minimises staff resources by fully inflating in only 2 minutes
- Hospital Grade Cover Thicker weight polyurethane fabric with increased bleach resistance to handle modern hospital infection control standards.
 Suitable for Sodium Hydrochloride solutions
- Adaptable Pump Pump can be switched between modes for typical users (up to 160 kg) and bariatric users (up to 450 kg) ensuring a single asset can meet the broadest spectrum of patient needs

Warranty Information

Cover: 2 Years

• Air Mattress: 2 Years

Pump: 2 Years



Compliance

IEC 60601-1-2

