

Aspire ALTITUDE Single Action Recliner



The Aspire Altitude Recliner is an electric recliner chair designed to provide comfortable seating for long periods of time. It is a truly mobile and flexible seating solution, with 75 mm castors facilitating smooth user transport plus a lithium-ion battery that can power the recliner all day when not within reach of power points.

- Use all day without plugging in the power with its inbuilt battery (requires recharging)
- Easy relocation of the chair between residents or facility locations with smooth-rolling 75 mm Tente castors (rear locking)
- Reduce shear and promote optimal skin microclimate with a multi-stretch vapour-permeable seating surface designed for users at risk of pressure injury
- Medical grade vinyl on the outer chair surface for improved longevity and infection control
- Removable armrests facilitate safe lateral transfers
- Enhanced comfort and immersive qualities with a three-tiered backrest cushion
- Optional oncology armrests (CHP227715) provide wide forearm support to enhance patient comfort and assist with effective administration of treatment such as blood transfusions
- Optional armrest table (CHP227720) allows for placement of personal items or medical instruments during treatment
- Some suitable applications include aged care facilities, rehabilitation facilities, hospitals, blood and treatment clinics, care in the home, and other independent in-home use



Contract approved

Support @ Home

Function: Sitting furniture

Category: Under advice

Click [here](#) for more information

Warranty Information

- Components: 2 Years*

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP227705	Aspire ALTITUDE Single Action Recliner - Black/Blue	800 mm	710 mm	1200 mm	150 kg	520 mm	510 mm	560 mm	560 mm	710 mm