

Rolyan SoftPro Functional Resting Hand



- Treats moderate contractures of the wrist/hand/fingers
- Easily adjustable for ulnar radial deviation
- Finger separator and ulnar drift strap available separately to allow custom fitting

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width
FBS301495LL	Rolyan SoftPro Functional Resting Hand - Left - Large	20.3 - 24.1	8.9 - 10.2
FBS301495LM	Rolyan SoftPro Functional Resting Hand - Left - Medium	15.2 - 20.3	7.6 - 8.9
FBS301495	Rolyan SoftPro Functional Resting Hand - Right - Small	Up to 15.2	Up to 7.6
FBS301495RM	Rolyan SoftPro Functional Resting Hand - Right - Medium	15.2 - 20.3	7.6 - 8.9
FBS301495RL	Rolyan SoftPro Functional Resting Hand - Right - Large	20.3 - 24.1	8.9 - 10.2
FBS301495LS	Rolyan SoftPro Functional Resting Hand - Left - Small	Up to 15.2	Up to 7.6