

Kerry Ezy Lift Patient Lifter



- Mobile patient lifter for transfer to a seated or lying position
- Ideal for use with weight scale or pivot frame
- Large load capacity of 225 kg
- Robust epoxy powder coated frame construction makes it ideal for aged care facilities
- Dual speed hand control operation
- Emergency lowering and raising
- Over-sized push handles for excellent manoeuvrability
- Electric leg gating
- Suits a large range of patient slings
- Supplied with battery and charger



Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 2 Years

| Code | Display Name | Total Length | Total Width | Total Height | Lift Range | Leg Gait Range | Electric Leg | Safe Working Load | Material |
|-----------|-------------------------------|--------------|-------------|----------------|---------------|----------------|--------------|-------------------|----------|
| LSS379800 | Kerry Ezy Lift Patient Lifter | 1440 mm | 690 mm | 1315 - 1730 mm | 485 - 1730 mm | 540 - 1000 mm | True | 225 kg | Steel |