

Kerry Ezy Lift Patient Lifter



- Mobile patient lifter for transfer to a seated or lying position
- Ideal for use with weight scale or pivot frame
- Large load capacity of 225 kg
- Robust epoxy powder coated frame construction makes it ideal for aged care facilities
- Dual speed hand control operation
- Emergency lowering and raising
- Over-sized push handles for excellent manoeuvrability
- Electric leg gating
- Suits a large range of patient slings
- Supplied with battery and charger



Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Electric Leg	Safe Working Load	Material
LSS379800	Kerry Ezy Lift Patient Lifter	1440 mm	690 mm	1315 - 1730 mm	485 - 1730 mm	540 - 1000 mm	True	225 kg	Steel