

## 2D Alice Bookshelf



- A modern bookcase with coastal-inspired styling featuring ample storage
- Features four spacious shelves for books and decorative items
- Generous cupboard space keeps clutter out of the way
- Light-tone teak finish that will seamlessly integrate into any living space
- Tall yet slim shape that will take up less floor space
- Dimensions: 1900mm (H) x 885mm (L) x 400mm (W)
- Pictured in Sublime Teak, colour/material to be advised at time of order



### Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height
FUL053110	2D Alice Bookshelf	885mm	400mm	1900mm