

## **Etac Hi-Loo Toilet Seat Raiser with Arm Support**



- The Etac Hi-Loo Toilet Seat Raiser with Fixed Mounting and Arm Supports is a practical two in one solution.
- Ergonomic smooth seat with recess at front and back for convenient intimate hygiene access.
- Securely mounted, it has two large friction elements for firm contact with the toilet.
- Detachable arm supports that swing up or can be permanently fixed in folded down position.

0			
(V)	Contracted	DVA	product

Ocontract approved

## Support @ Home

**Function:** Toileting **Category:** Prescribed

Click **here** for more information

Code	Display Name	Total Height	Safe Working Load
BTT162005	ETAC Hi-Loo 60mm	60 mm	150 kg
BTT162010	ETAC Hi-Loo 100mm	100 mm	150 kg
BTT162015	ETAC Angled Seat Hi-Loo	100 - 60 mm Angled	150 kg