

## Pride C101 Riser Recliner Powerlift - Single Motor



Pride's NEW Riser Recliner delivers superb value, excellent quality and luxurious comfort. The Riser Recliner features quiet and smooth operation and stylish design to add to any home decor.

- Rise & Recline (single motor - foot raise with recline function)
- Medium to large lift chair
- Engineered hardwood frame
- Emergency battery backup system in case of power failure

### Warranty Information

- 2 Years

**PRIDE**  
MOBILITY  
LIVE YOUR BEST®

✓ Contracted DVA product



**ndis**

**Item Code:**

05\_180903180\_0103\_1\_2

Click [here](#) for more information

### Compliance

ISO13934-1 / ISO13937-1 / BS6807  
/ BS7175

Code	Display Name	RRP	Box Quantity	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	User Height	Safe Working Load	Colour
CHP197050	C101 Riser Recliner - Chocolate	1719	1	510 mm	510 mm	700 mm	470 mm	1630 - 1830 mm	147 kg	Pecan Fabric