

## Liko MasterVest Walking Sling Leg Harness Mod 66 - Poly



- Leg Harness for MasterVest [Liko Mastervest](#)
- Versatile training aid that facilitates toilet visits, standing training, gait training, mobility training and placing patients into standing positions
- Reduces strain on caregivers while helping them deliver physical training
- Transfers to an from the toilet can be accomplished easily and smoothly, reducing strain
- Lifts around the upper body and partly under the arms, and is best suited for patients that can support some of their body weight
- Available in Polyester



### ndis

**Item Code:**

05\_123603166\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

### Warranty Information

- 1 Year

### Compliance

AS/NZS ISO 10535:2011



Code	Display Name	Safe Working Load	Material
LSS394310	MasterVest Walking Sling Harness Mod 66 - Polyester - XSMALL / SMALL	200 kg	Leg Harness (Optional)
LSS394320	MasterVest Walking Sling Harness Mod 66 - Polyester - MEDIUM / LARGE	200 kg	Leg Harness (Optional)
LSS394330	MasterVest Walking Sling Leg Harness Mod 66 - Polyester - XLARGE / XXLARGE	200 kg	Leg Harness (Optional)