

Kerry Access Toileting Sling - Polyester



- Attaches to hoist to enable a carer to lift and manoeuvre a less mobile user for toileting
- User must have good trunk control and be partial weight bearing
- Quick and easy to fit
- Colour coded attachment and positioning straps
- Note: should not be used for non-weight bearing patients
- Machine washable
- FIM Scores: 2 and 3



ndis

Item Code:

05_123603163_0105_1_2

Sling – Composite:

05_421_0105_1_2

Click [here](#) for more information

Warranty Information

- 1 Year

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS372500	Access Toileting Sling - Small	250 kg	Polyester
LSS372600	Access Toileting Sling - Medium	250 kg	Polyester
LSS372700	Access Toileting Sling - Large	250 kg	Polyester
LSS372710	Access Toileting Sling - Extra Large	250 kg	Polyester