

Queens Cutlery



 A wide range of stainless steel utensils mounted into a built-up plastic handle to help those with weak or restricted grip, the handle is cylindrical to assist with gripping

• The items are hygienically sealed, easy to clean and are dishwasher safe

Handle Length: 98 mmHandle Diameter: 30 mmProduct Weight: 45 g

Warranty Information

• 1 Year

Code	Display Name
DLK275760	Knife
DLK275770	Fork
DLK275780	Dessert Spoon
DLK275790	Junior Spoon

DLK275870 Cutlery Set

Ontracted DVA product

Ocontract approved

ndis

Item Code:

05_150900319_0103_1_2 Click **here** for more information

Compliance

No applicable compliance standard