

## Queens Cutlery



- A wide range of stainless steel utensils mounted into a built-up plastic handle to help those with weak or restricted grip, the handle is cylindrical to assist with gripping
- The items are hygienically sealed, easy to clean and are dishwasher safe
- Handle Length: 98 mm
- Handle Diameter: 30 mm
- Product Weight: 45 g

### Warranty Information

- 1 Year

- ✔ Contracted DVA product
- ✔ Contract approved

### Support @ Home

**Function:** Eating and drinking

**Category:** Low risk

Click [here](#) for more information

### NDIS

**Item Code:**

05\_150900319\_0103\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	Unit Weight
DLK275760	Knife	45 g
DLK275770	Fork	45 g
DLK275780	Dessert Spoon	45 g

Code	Display Name	Unit Weight
DLK275790	Junior Spoon	45 g
DLK275870	Cutlery Set	-