

Denver Chair



- Designed to bring pressure management, posture support and comfort
- Dual motors provide leg elevation and tilt in space positioning as well as adjustable recline angle. The tilt in space combined with back recline can reduce sliding from the chair, while the riser function helps the user to a standing position to promote easier independent transfers
- Full pressure care materials throughout the sitting surface, including a high grade memory foam seat cushion and waterfall back rest. All areas of the Denver are able to be cleaned and the cushion can be removed
- Adjustable in seat depth, back rest configuration and has the option of lateral supports to help maintain posture
- For multi-user environments, hospitals and equipment loan stores, the Denver is completely modular and all parts can be replaced

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Sitting furniture

Category: Under advice

Click [here](#) for more information

Warranty Information

- 7 Years Frame Warranty
- 1 Year Motor Warranty

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP197832	Chair Pressure Relief - Denver Recliner Chair - 450mm - Electric	1700 mm	700 mm	1090 mm	137 kg	457 - 559 mm	450 mm	400 - 500 mm	400 - 500 mm	690 mm