

iCare - Mattress Cover



- Fully wipeable with welded seams and drip cover over zipping
- Perfect for at home or in facilities where infection control is required
- Easy to wipe down and clean
- Long Single: Fits mattress size 900 x 2030 mm
- King Single: Fits mattress size 1030 x 2030 mm

 Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Total Width
BEM469710	Mattress Cover - Long Single	2030 mm	900 mm
BEM469715	Mattress Cover - King Single	2030 mm	1070 mm