

Stethoscope Luxury Dual Head



- Well designed economical stethoscope
- Lightweight aluminum reversible head & chest piece
- Non chill ring for patient comfort

Warranty Information

- 1 Year

✔ Contract approved

Support @ Home

Function: Not on S@H

Category: Not on S@H

Click [here](#) for more information

Code	Display Name
GMD300100	Stethoscope Luxury Dual Head - Black
GMD300110	Stethoscope Luxury Dual Head - Magenta
GMD300120	Stethoscope Luxury Dual Head - Navy Blue