

Raised Toilet Seat with Swing-Back Arms - Wide



- Width Between Arms: 430 mm
- Unique design to ease the use of a toilet by heightening the seat and adding arm rests to hold
- Fits most Toilets
- Easy to Install
- Easy to Clean
- Arm Swings Back for Convenience

✓ Contract approved

Support @ Home

Function: Toileting

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BTT150460	Raised Toilet Seat with Swing-Back Arms - Wide	490 mm	530 mm	450 - 550 mm	120 kg