

Oppo 1032 Knee Support Post Operative



- Reduces swelling and pain following arthroscopic surgery and enhances protection of periarticular structures
- Provides comfortable and controlled compression with adjustable hook and loop closures
- Front opening design for easy fit and comfort

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Knee Circumference
FBS301634L	Oppo - Knee Hinged - Thermal Wrap Post Op - O/Patella - Lge	406 - 451 mm
FBS301634XL	Oppo - Knee Hinged - Thermal Wrap Post Op - O/Patella - XL	451 - 502 mm
FBS301634M	Oppo - Knee Hinged - Thermal Wrap Post Op - O/Patella - Med	375 - 406 mm
FBS301634S	Oppo - Knee Hinged - Thermal Wrap Post Op - O/Patella - Sml	375 - 406 mm