

Regent Bedside Cabinet



- Classic design adds an elegant touch to any bedroom
- A versatile piece that can be used in bedroom, living area or office
- Functional and slimline design perfect for smaller spaces
- Features four drawers with ample storage for personal items and documents
- Lockable drawer can be used to store medications and valuables
- Ample space on top for a lamp, phone and book
- Dimensions: 705 mm (H) x 475 mm (L) x 450 mm (W)
- Pictured in Siam Teak, colour/material to be advised at time of order



Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height
BEA828430	Regent Bedside Cabinet	475mm	450mm	705mm