

Aidacare Paediatric Cart



- The Aidacare Paediatric Crash Cart has sleek lines, and strong construction
- Aidacare carts are designed to offer an attractive, practical and durable cart with multiple configurations to enable the efficient organisation and delivery of medication and equipment
- Robust panels to minimise dents and scratches
- Large work area with slide out work surface are standard
- Easy roll 125 mm castors
- Designed to accommodate the Broselow System of colour coded drawers.
- Ensures doctors and technicians can quickly identify critical supplies in emergencies – when records matter
- Fitted standard with a keyless break-away lock and comprehensive range of standard accessories
- The simple breakaway lock offers easy access, combined with essential security
- **Standard Inclusions:** Defibrillator Shelf, Oxygen Bottle Holder, Extended Sliding Surface, Waste Bin
- Customisation to suit a range of applications are also available upon request
- Comprehensive range of standard accessories available
- Standard drawer configuration Paediatric Crash Cart:
 - 8 small 3' drawers (7 cm)
 - 1 medium 6' drawer (16 cm)



Accessories

Click on the cart icon to display the range of optional accessories

Warranty Information

- Structural: 5 Years Warranty
- Accessories: 1 Year Warranty

Code	Display Name	Total Length	Total Width	Total Height	Colour
WAT909690	Aidacare Paediatric Cart - Large	790 mm	610 mm	1100 mm	Multi