

Seville Floorline Folding Bed



- Offers the user and carer an easy to install and transportable floor line bed which is ideal for use in most care environments
- Integrated 200 mm bed extension
- Lock out handset
- Battery backup
- Trendelenburg (15°)
- Configurations: 4 Function - Electric high/low, backrest, knee bend and Trendelenburg
- Recommended mattresses: Lifecomfort Deluxe AC BEM046510



Compliance

No applicable compliance standard

Warranty Information

- 10 Years Frame
- 2 Years Electrics

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Number of Functions
BEB045335	Seville Floorline Folding Bed	2340 - 2540 mm	910 mm	130 - 725 mm	225 kg	4