

## TributeWrap Wrist to Axilla



- TributeWrap is an adjustable compression garment with foam technology that nestles into the tissue to promote skin health. It supports easy self-management of lymphoedema and oedema during low activity by comfortably adjusting to your patient's unique shape and lifestyle
- TributeWrap helps maintain gains made during therapy while allowing the flexibility to support further reduction. Simply slide it on, adjust the straps, and cover with Sleep Sleeve for comfortable evening compression

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Supporting and stimulating blood circulation

**Category:** Low risk

Click [here](#) for more information

Code	Display Name
FBS310140	TributeWrap Wrist to Axilla Small Right
FBS310130	TributeWrap Wrist to Axilla Small Left
FBS310120	TributeWrap Wrist to Axilla Medium Right
FBS310110	TributeWrap Wrist to Axilla Medium Left
FBS310100	TributeWrap Wrist to Axilla Large Right
FBS310090	TributeWrap Wrist to Axilla Large Left