

Aspire ComfiMotion Activ Care Bed Rails



EXTENDABLE BED RAIL BEA830205



- Improves bed edge awareness and offers multiple grip points for repositioning and bed mobility
- Designed to be used in conjunction with Comfimation Activ Care Bed
- Extendable Version (BEA830205): Ability to extend to prevent falls and contract to allow safe bed egress and ingress
- Additional specifications
 - Height of rail from mattress deck: 565 mm
 - BEA830205 width of rail: 680 - 1060 mm
 - BEA830212 width of rail: 685 mm

 Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Code	Display Name
BEA830205	ComfiMotion Activ Care Bed - Extendable Bed Rail
BEA830212	ComfiMotion Activ Care Bed - Bed Rail