

Aidacare Aspire Inline Weigh Scale



- Designed to be fitted on Aidacare Aspire Patient Lifters
- Allows for an accurate weight measurement to be taken by a carer when a patient is suspended in a sling
- Large weight capacity of up to 320kg
- Easy to read LCD display

✓ Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load
LSS390680	Aspire In-Line Weigh Scale	320 kg