

Cura Alarms - Patient Activated Monitors



cura¹

- Patient activated alert monitor assists with falls prevention
- Attached to patients clothing, patient activates by pulling string which detaches magnet and sets off alarm
- Stand alone system is simple to use
- Durable and reliable
- Portable and lightweight
- Supplied with 9 volt battery

Warranty Information

- 1 Year

 Contract approved

ndis

Item Code:

05_222704219_0123_1_2

Personal Alarm – Composite:

05_427_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name
BEA007460	Pull String Monitor