

Lifting Hoist Heavy Duty (200kg) - Includes 1 Sling



A robust workhorse, the Aidacare Aspire 200 is designed to meet the rigorous demands of professional and facility environments

- Ergonomic Handset. The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- Powered Leg Function. Opens and closes the legs for lifting stability or accessibility
- Quick Release Yoke Attachment. One touch release and attach mechanism, designed for the busy caregiver. Padded yoke for added patient safety
- Sturdy steel frame with industry leading Linak electronics and warranty provision
- Safety Critical features include emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart
- Intelligent Diagnostics. State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- Quality Locking Castors. 2 x 100mm rear locking castors for safety and manoeuvrability

*** Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

Code

Display Name

H_LSS390620 Lifting Hoist Heavy Duty (200kg) - Includes 1 Sling

