

Aspire ComfiMotion Pocket Spring Mattress



- The ComfiMotion Pocket Spring mattress offers zoned support springs that respond and contour to the individual body shape to improve comfort, maximising support and sleep quality
- Transitional layered pocket springs arranged in weight activated support zones
- Integrated side-walls create roll resistance at the bed edge and provide extra rigidity to ease bed entry and exit
- Luxuriously soft pillowtop layer provides a plush immersive feeling complimented by the responsive pocket springs

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 10 Year Limited Warranty - Memory Foam
- 5 Year Limited Warranty - Pocket Spring

Compliance

BS7177:2008



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEM469800	Aspire ComfiMotion Pocket Spring Mattress - Portobello Grey (sides) - Long Single	2030 mm	920 mm	240 mm	180 kg
BEM469810	Aspire ComfiMotion Pocket Spring Mattress - Portobello Grey (sides) - King Single	2030 mm	1070 mm	240 mm	180 kg
BEM469820	Aspire ComfiMotion Pocket Spring Mattress - Portobello Grey (sides) - Long Double	2030 mm	1370 mm	240 mm	360 kg
BEM469830	Aspire ComfiMotion Pocket Spring Mattress - Portobello Grey (sides) - Queen	2030 mm	1520 mm	240 mm	360 kg