

Aspire Active Air Acute 8 Mattress Replacement



- The Active Air Acute 8 is specifically designed to meet the demands of Critical Care environments, maximising patient comfort and therapeutic outcomes
- 1 in 3 Cell Alternation Cycle with 24 Narrow Profile cells – Improves patient comfort by reducing space between cells and increasing patient support surface area
- Focused Heel Care – Removable heel cells facilitate rapid response to changes in heel therapy needs
- Dual Compressor Pump – Reduces noise and efficiently manages patient therapy at highest patient weight loads (250 kg)
- Hospital Grade Cover – Thicker weight polyurethane fabric with increased bleach resistance to handle modern hospital infection control standards. Suitable for Sodium Hydrochloride solutions
- Lateral turning through integrated air cells, allowing the mattress to turn a patient without the need to add an after-market turning device
- Drastically reduces the moving and handling stresses on Care Staff by rotating the patient up to 40 degrees. This technology is ideal for:
 - Alleviating asymmetric or laterally located pressure areas
 - Simplified in-bed bathing
 - Rapid changing of linen
 - Perfect positioning of manual handling slings and slide sheets

✓ Contract approved

Compliance

IEC 60601-1-2

Warranty Information

- Cover: 2 Years
- Air Mattress: 2 Years
- Pump: 2 Years

| Code | Display Name | RRP | Box Quantity | Total Length | Total Width | Total Height | Safe Working Load | Pressure Care Rating |
|------------|--|------|--------------|--------------|-------------|--------------|-------------------|----------------------|
| PTM639400T | Active Air ACUTE 8 w/ turnassist – Single | 4919 | 1 | 2000 mm | 900 mm | 200 mm | 30 kg - 250 kg | Very High Risk |
| PTM639410T | Active Air ACUTE 8 w/ turnassist – King Single | 6319 | 1 | 2000 mm | 1050 mm | 200 mm | 30 kg - 250 kg | Very High Risk |