

Teal Roku Bed with Side Storage



A secure, extremely durable one-piece rotationally moulded bed for challenging environments.

- Smooth, high impact polyethylene, continuously moulded for maximum durability
- Underbed storage area, ideal for client personal effects
- Surface grooves facilitate easier cleaning and improved infection control
- Seamless and stain-resistant design improves longevity
- Integrated mattress retainers prevent mattress movement
- Smooth, round profile edge ensures enhanced safety and security for users and carers
- 100% recyclable material
- Pairs with the Roku® range of bedroom solutions
- Tested and compliant to BS EN 597-1, BS EN 597-2, BS6807:2006 clause 9 & BS 7177:2008

✓ Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
BEB050520	Teal Roku Bed w/side Storage - (Choose colour)	2120 mm	1005 mm	510 mm	320 kg	71 kg