

BodyAssist Elastic Strap with Loop



- Multi purpose elastic strap with a loop sewn at one end for effective wrapping
- Ideal for wrists, arms, ankles, calves and forearms
- Designed to stabilise and assist healing of injured soft tissue
- One size fit



Warranty Information

- 1 Year

| Code | Display Name | Total Length | Colour |
|-----------|-------------------------------------|--------------|--------|
| FBS283540 | Elastic Strapping - w/ Loop - Black | 700 mm | Black |
| FBS283550 | Elastic Strapping - w/ Loop - Beige | 700 mm | Beige |
| FBS283555 | Elastic Strapping - w/ Loop - White | 700 mm | White |