

## Mobile IV Stand



- 2 point “pig tail” IV bag attachment
- Adjustable pole with pneumatic dampener mechanism to retard adjustment speed
- Stainless steel 25 mm diameter upright can accept 1 x syringe pump
- Stable 5 star stainless steel base 500mm diameter

 Contract approved

### Compliance

No applicable compliance standard

### Warranty Information

- 5 Years

| Code      | Display Name  | Total Length | Total Width | Total Height      | Unit Weight |
|-----------|---|--------------|-------------|-------------------|-------------|
| WAT820130 | IV pole - Mobile - 100mm Grey Rubber Castors                | -            | -           | 1350 - 2250 mm    | -           |
| WAT820140 | IV pole - Mobile - 100mm Grey Castors + Push Ring + 4 Hooks | -            | -           | 1350 - 2250 mm    | -           |
| WAT820170 | Push Handle Ring - To Suit IV Pole & Syringe Pump Range     | -            | -           | -                 | -           |
| WAT820145 | IV Pole - 4 Infusion with 6x3 Litre bags                    | 600 mm       | 600 mm      | 1550 mm - 1950 mm | 80 kg       |