

TED - Anti-Embolic - Below Knee Stockings



- TED knee high stockings are designed to reduce venous stasis below the knee when thigh length styles are medically contraindicated
- Knee length with soft top band to minimise slippage without constricting circulation
- Available in 6 sizes, with choice of either regular or long fit

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name	Colour
PTS641570	TED Anti-embolic Stockings - Below Knee - Toe Hole - Sml / Reg	White
PTS641580	TED Anti-embolic Stockings - Below Knee - Toe Hole - Med / Reg	White
PTS641593	TED Anti-embolic Stockings - Below Knee - Toe Hole - XXXL / Reg	White
PTS641600	TED Anti-embolic Stockings - Below-Knee - Toe Hole - Sml / Long	White
PTS641610	TED Anti-embolic Stockings - Below-Knee - Toe Hole - Med / Long	White
PTS641620	TED Anti-embolic Stockings - Below-Knee - Toe Hole - Lge / Long	White
PTS641621	TED Anti-embolic Stockings - Below-Knee - Toe Hole - XL / Long	White

Code	Display Name	Colour
PTS641622	TED Anti-embolic Stockings - Below-Knee - Toe Hole - XXL / Long	White
PTS641623	TED Anti-embolic Stockings - Below-Knee - Toe Hole - XXXL / Long	White