

## 200S Sling Bar - Clips



- Designed to convert the Aidacare Aspire 200S Standing Lifter from a loop sling attachment to a clip sling attachment
- Facilitates the use of Aspire Standing Slings with Clips
- Easily retrofittable to new and existing 200S lifters

### Warranty Information

- 2 Years Frame

✔ Contract approved

### Support @ Home

**Function:** Lifting persons

**Category:** Prescribed

Click [here](#) for more information



Code	Display Name	Safe Working Load
LSS390655	200S Sling Bar with Clips	200 kg