

Aspire Low Friction Sling with Head Support - Parasilk - Black Clips



- Clip fitting option to fit sling to a pivot frame
- Australian standards tested and registered with the TGA
- Parasilk material, with low friction properties, for maximum patient comfort during transfers and easy positioning
- Quality - manufactured in the UK using European materials, computer guided stitching and CNC fabric cutting
- Durable keyhole clip attachment points for safety and security during transfers
- Serial Numbered - slings are embroidered with a unique serial number to maintain traceability and compliance at all times
- Service & Maintenance - service record labels ensure you can quickly and easily identify service history
- Easy Sling Identification - all slings are colour coded by size and have reinforced edging for longevity
- Easy Clean - machine washable up to 85°C temperature
- Parasilk leg strap lining for easy positioning
- For FIM score 1 or 2



Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 1 Year

Code	Display Name	RRP	Box Quantity	Safe Working Load	Material
LSS390960B	Deluxe Low Friction Sling - Parasilk - Black Clips - Small	559	1	200 kg	Parasilk

Code	Display Name	RRP	Box Quantity	Safe Working Load	Material
LSS390965B	Deluxe Low Friction Sling - Parasilk - Black Clips - Medium	559	1	200 kg	Parasilk
LSS390970B	Deluxe Low Friction Sling - Parasilk - Black Clips - Large	559	1	300 kg	Parasilk
LSS390975B	Deluxe Low Friction Sling - Parasilk - Black Clips - Extra Large	559	1	300 kg	Parasilk