

Seat Walker Classic (130kg)



- Ultra-soft hand brakes with palm-ball lock effortlessly, offering increased comfort for users with arthritis or wrist/hand joint pain
- Handle height adjusts quickly with hand wheels for optimal sizing
- Rustproof, robust aluminium frame provides durability and stability and whilst remaining very light to manoeuvre and transport
- Easy folding mechanism facilitates storage and transportation
- Anatomically designed, latex-free handles allow comfortable grip especially for users with arthritis
- Curved and padded backrest offers comfortable seated support and flips up/removes easily for folding/storage
- Padded seat flips up easily for access to storage bag (supplied) or for folding
- Swiveling front castor wheels increase maneuverability
- 8 inch diameter solid, puncture-proof wheels provide increased maneuverability over various surfaces
- 10 Easy wipe-clean surfaces



*** Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

| Code | Display Name |
|---------------|-----------------------------|
| H_WAF705700BL | Seat Walker Classic (130kg) |