

Linen Bag - Drawstring



- Laundry bags in various colours with drawstring option to keep linen secure when transporting from patients rooms to laundry facilities
- Heavy duty draw cord
- Industrial strength cord lock
- Robust bottom handle
- Fits most trolleys

✔ Contract approved

Compliance

No applicable compliance standard

Warranty Information

- 1 Year

Code	Display Name	Material
LAC353130	Drawstring Linen Bag - White	Polyester
LAC353140	Drawstring Linen Bag - Blue	Polyester
LAC353160	Drawstring Linen Bag - Green	Polyester
LAC353170	Drawstring Linen Bag - Red	Polyester
LAC353180	Drawstring Linen Bag - Yellow	Polyester
LAC353190	Drawstring Linen Bag - Purple	Polyester
LAC353150	Drawstring Linen Bag - Brown	Polyester

Code	Display Name	Material
-------------	---------------------	-----------------