

## Aspire COVE Pressure Relief Chair



- The Aspire Cove Chair is a mobile chair designed to improve postural, pressure care and comfort outcomes for users requiring medium to high level support. Users can be easily transported in comfort to facilitate inclusiveness and social interaction
- **Adjustable Pressure Seat Cushion:** A hybrid air-foam seat cell offers convenient pressure adjustment to carers and maintenance staff. Reinflation is as simple as opening and closing a valve
- **Market Leading Upholstery:** The patient surface is genuine DARTEX made in the UK. This breathable highly elastic surface maximises skin microclimate and pressure redistribution. External surfaces are medical grade coated fabrics from Warwick Australia. These hardy yet stylish patterns offer an appealing aesthetic and blend seamlessly in with modern Aged Care interior designs
- **Depth Adjustable Seat:** Backrest can be adjusted forward/back to achieve the correct seat depth and fit for each user, also maximising lower leg support when elevated
- **Infinite Tilt and Leg Positioning:** Carers can easily adjust the Tilt position through a 40° range to effectively redistribute sitting pressure away from vulnerable bony prominences. Leg Raise operates independently of Tilt to allow infinite positioning combinations and maximise lower leg support. Both functions are gas assisted for the ease of carers
- **Infection Control:** Medical grade upholstery with removable seating elements simplify infection control processes. The biostatic patient surface prevents bacteria and virus ingress
- **Contoured Head Support Pillow:** Comfortable head and neck support adjusts easily to achieve the desired position for each user.
- **Transitional Support Armrests:** The immersive padded armrests track arm positioning throughout the tilting range. This promotes user relaxation and

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Sitting furniture

**Category:** Under advice

Click [here](#) for more information

### Accessories

Click on the cart icon to display the range of optional accessories

confidence when being repositioned

- **Ergonomic Carer Push Handles:** Multiple grip positions offer improved carer ergonomics throughout the tilting range, making for easy operation and user repositioning
- **Anatomical Legrest Movement:** As the legrest raises, it extends outwards tracking knee joint range of motion and maximising lower limb support
- **Tool Free Adjustment:** Critical adjustments can be made easily without the need for tools. This improves speed of care delivery and maximises the use of seating assets in multi-user environments
- **Ease of Movement:** High performance TENTE Levina castors reduce push/pull forces during transportation. This improves workplace health and safety outcomes for care staff and family members
- **Complex Pressure Care Compatibility:** The layered seat structure accommodates after-market pressure care cushions such as ROHO replacing the air-foam cushion where necessary, without compromising chair performance or infection control outcomes
- **Advanced Immersion Backrest:** This backrest combines contouring with a plush and immersive surface to provide greater lateral support when upright sitting and caters more directly for users with kyphosis ([optional accessory available here](#))
- **Utility Table:** Offers a convenient surface for refreshments, books, electronic devices and the like. Rotates out of the way for stowage or for clear seat entry and exit. Mounted to the non-moving armrest frame, the tray stays horizontal while a user is tilted and repositioned ([optional accessory available here](#))
- **Weight Bearing Footrest:** The fully retractable footrest protects the user's feet during transportation, redistributes pressure away from bony prominences during prolonged upright sitting and serves as a step to assist with seat entry and exit ([optional accessory available here](#))

## Warranty Information

- Fabric: 1 year
- Frame: 3 Years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP198960	Aspire COVE Pressure Relief Chair - Medium - Grey/Andreas Pumice	860 - 1640 mm	780 mm	1270 mm	160 kg	Grey Andreas Pumice	500 - 580 mm	550 mm	540 mm	470 mm	740 mm
CHP198965	Aspire COVE Pressure Relief Chair - Medium - Navy/Bronx Steel	860 - 1640 mm	780 mm	1270 mm	160 kg	Navy Bronx Steel	500 - 580 mm	550 mm	540 mm	470 mm	740 mm
CHP198970	Aspire COVE Pressure Relief Chair - Large - Grey/Andreas Pumice	860 - 1690 mm	830 mm	1300 mm	160 kg	Grey Andreas Pumice	500 - 580 mm	600 mm	540 mm	470 mm	740 mm
CHP198975	Aspire COVE Pressure Relief Chair - Large - Navy/Bronx Steel	860 - 1690 mm	830 mm	1300 mm	160 kg	Navy Bronx Steel	500 - 580 mm	600 mm	540 mm	470 mm	740 mm
CHP198978	Aspire Cove Pressure Relief Chair - Medium - Teal Vinyl	1640 mm	780 mm	1270 mm	160 kg	Teal Vinyl	500 - 580 mm	550 mm	540 mm	540 mm	740 mm
CHP198979	Aspire Cove Pressure Relief Chair - Large - Teal Vinyl	1690 mm	830 mm	1300 mm	160 kg	Teal Vinyl	500 - 580 mm	600 mm	540 mm	540 mm	740 mm