

Careline Huntington Bed



- Made in Germany, the Careline Huntingtons bed is of the highest quality. The beds are sturdy in construction and built to resist constant movements by the user. The mattress is also waterproof and designed to reduce the risk of shearing
- The bed offers high low functions, backrest and leg support in addition to being fully padded for the protection of the user. This padding extends to the head and footboards and also the rails are padded
- The padded rails are electrically operated to minimize risk to users and also carers. The castors are lockable by a foot control
- Bed Height Range: 410 - 610 mm
- Backrest incline/decline: 0 - 55 degree
- Knee bend incline/decline: 0 - 40 degree
- Mattress Dimensions: 2000 x 900 mm

Warranty Information

- 1 Year



Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_181210197_0103_1_2

Power Bed – Composite:

05_423_0103_1_2

Click [here](#) for more information

Compliance

IEC 60601-2-52:2009

Code	Display Name	Total Width	Safe Working Load
BEB046010	Huntington's Bed - Electric	900 - 1100 mm	150 kg