

Pride C1 Powerlift Recliner - Single Motor



PRIDE
MOBILITY
LIVE YOUR BEST®

- Rise & Recline (single motor - foot raise only with recline function)
- Petite to small lift chair with fully removable upholstery components
- Engineered furniture hardwood frame
- Emergency battery backup system in case of power failure
- Common Colours: Arctic Blue and Chocolate. Also available in Oatmeal (with upcharge)
- For DVA: Colour subject to availability

Warranty Information

- 2 Years

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_180903180_0103_1_2

Click [here](#) for more information

Compliance

ISO13934-1 / ISO13937-1 / BS6807
/ BS7175

Code	Display Name	RRP	Box Quantity	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	User Height	Safe Working Load	Colour
CHP197209	C1 Recliner - Arctic Blue	1949	1	510 mm	460 mm	660 mm	410 - 460 mm	< 1630 mm	113 kg	Arctic Blue Fabric
CHP197211	C1 Recliner - Chocolate	1949	1	510 mm	460 mm	660 mm	410 - 460 mm	< 1630 mm	113 kg	Chocolate Fabric