

Teal Myla Security Bed

A secure freestanding bed designed for challenging environments. Smooth, round profile edges enhances safety and security.



- Durable MDF vinyl wrapped soft edge top (25mm thick)
- Channeled platform allows for air to circulate under mattress
- Waterfall edging around the entire perimeter prevents fluid ingress and facilitates easier cleaning for improved infection control
- Sloping sides enable easy staff access for cleaning or assisting a user to aid sit-to-stand movements
- Integrated mattress retainers prevent mattress movement
- Pair your Myla with a mattress solution, contact Aidacare to customise your solution
- Tested and compliant to BS EN 597-1, BS EN 597-2, BS6807:2006 clause 9 & BS 7177:2008

 Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
BEB050500	Teal Myla Security Bed - (Choose colour)	1950 mm	1020 mm	385 mm	320 kg	78 kg