

MOTomed Loop Light.LA



RECK
MOTomed®

The MOTomed Loop brings movement therapy into the future. Through the new design and its digital intelligence the user can experience passive or active training in a completely new way. With the new and innovative height adjustment, users can adjust the position of the MOTomed loop without using tools. The Light model provides high quality primary care with all the common functions required from a movement therapy trainer - including Spasm Detection program for maximum user safety.



- The LA comes with both upper limb and lower limb trainers installed
- Soft Grip safety foot shells with fastening straps (pair, incl. pedal cranks)
- Soft Grip Handles (pair, incl. pedal cranks)
- Gas spring support for height adjustment
- Special stand with reduced stand height
- Leg guides with plastic-coated calf shells (pair)
- QuickFix foot fastening system (pair)
- Hand fixation by wrist cuff (piece)

| Code | Display Name | Total Length | Total Width | Total Height | Safe Working Load | Unit Weight |
|-----------|--|--------------|-------------|--------------|-------------------|-------------|
| FPM101140 | MOTomed Loop.LA Light - Serial Interface, Gas Spring, Stand, Leg guides, Foot & Hand Fix | 70 cm | 60 cm | 107 - 122 cm | 135 kg | 36 kg |
| FPM101141 | MOTomed Loop LA Light Basic - Leg Guides | - | - | - | - | - |