

Slingbar Mini 220



- SlingBar Mini 220 is an especially small sling bar which has been developed mainly for patients who prefer to sit 'snugly bound' in their sling
- Normally, the explanation is that the sling should be as tight as possible, so that it feels extra safe and secure

Warranty Information

- 1 year

 Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011



Code	Display Name	Safe Working Load
LSS393460	Spreader Bar - Liko SlingBar Mini 220	205kg