

Rehband 5055 Acro Comfort Shoulder Support



- Stabilises the shoulder joint
- Stretch flex system is anatomically designed for optimal function
- Excellent comfort through the temperature regulating features of comfortemp®
- Fits left or right

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Under Arm/Scapula Circumference
FBS301257	Rehband - Shoulder Support - Acro Comfort - X Sml	240 - 280 mm
FBS301257L	Rehband - Shoulder Support - Acro Comfort - Lge	350 - 380 mm
FBS301257M	Rehband - Shoulder Support - Acro Comfort - Med	320 - 350 mm
FBS301257S	Rehband - Shoulder Support - Acro Comfort - Sml	280 - 320 mm
FBS301257XL	Rehband - Shoulder Support - Acro Comfort - XL	380 - 420 mm