


Lifecomfort Bariatric Mattress (300kg)



- Pressure reducing bariatric mattress
- Extra weight capacity up to 300 kg
- Four layer foam support 200 mm high for bariatric use with maximum comfort and support
- Full size memory foam overlay
- Foam side walls for extra stability on the edge of the mattress
- Two way waterproof stretch cover
- PVC bottom cover
- Four sided zip for easy cover removal
- Stitched zipper cover with waterproof flap for excellent hygiene
- Independently tested and certified for durability and weight rating
- PUPPS 2 (2004) compliant



*** Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

| Code | Display Name |
|-------------|--|
| H_BEM046700 | Lifecomfort Bariatric Mattress (300kg) |