

HASOMED RehaMove2 Accessories



- A group of accessories and replacement parts for the RehaMove2.



HASOMED+

Code	Display Name	Total Length	Total Width
FPM100110	RehaMove2 - Electrodes - Rectangle	90 mm	50 mm
FPM100120	RehaMove2 - Electrodes - Large Rectangle	130 mm	75 mm
FPM100130	RehaMove2 - Electrodes - Oval	64 mm	40 mm
FPM100175	RehaMove2 - Electrode 8 Channel Replacement Lead (Right) -	-	-
FPM100180	RehaMove2 - Electrode 8 Channel Replacement Lead (Left) -	-	-
FPM100190	RehaMove2 - Power Supply	-	-
FPM100200	RehaMove2 - Carry bag	-	-

Code	Display Name	Total Length Total Width
------	--------------	--------------------------